


# DAD Detector Model LaChrom Elite L-2450

## FLOW CELL, 10 MM



Durchflußzelle Standard  
Bestell Nr. 890-2394

☒ L-2450

## TUBE INLET ASSY for standard flow cell



Einlass Tube Standard  
Bestell Nr. 890-2288

☒ L-2450

## TUBE OUTLET ASSY for standard flow cell



Auslass Tube Standard  
Bestell Nr. 890-2289

☒ L-2450

## LAMP D2 (Original Hitachi)



Original D2 Lampe  
Bestell Nr. 890-2430

☒ L-2450

## Lamp D2 (Hamamatsu 1000h)



Hamamatsu D2 Lampe  
Bestell Nr. L-5257

☒ L-2450

## L2D2 (Hamamatsu 2000h)



L2D2 2000h Long-Life  
Bestell Nr. L-6901

☒ L-2450

## LAMP Wolfram



Wolfram Lampe  
Bestell Nr. 890-2527

☒ L-2450

## CELL PACKING)



Distanz Ring  
Bestell Nr. 890-2166

☒ L-2450

## SPACER for Flow Cell



Distanz Ring  
Bestell Nr. 060-3076

☒ L-2450

## O-RING for Flow Cell



O-Ring  
Bestell Nr. L456005

☒ L-2450