

### STATOR FOR 7010/7120/7125 RHEODYNE



Stator Injektion Valve  
Bestell Nr. 7010-040

L-7200  L-7250

### ROTOR FOR RHEODYNE VESPEL



Rotor Vespel -pH 1-8  
Bestell Nr. 7010-039

L-7200  L-7250

### ROTOR FOR RHEODYNE TEFZEL



Rotor Tefzel - pH 1-14  
Bestell Nr. 7010-071

L-7200  L-7250


### PIPE Capillary Port-Injector



Injektions Port  
Bestell Nr. 810-3547

L-7200  L-7250

### PIPE DRAIN Valve-Wasteport



Waste Port  
Bestell Nr. 810-3548

L-7200  L-7250

### SEAL for Injection-Port



Dichtung Inj.port  
Bestell Nr. 080-3040

L-7200  L-7250

### STATOR for Dilutor Valve



Stator Spülventil  
Bestell Nr. 810-3097

L-7200  L-7250

### DILUTOR VALVE SEAL



Rotor Spülventil  
Bestell Nr. 810-3095

L-7200  L-7250

### BEARING



Bearing Spülventil  
Bestell Nr. 810-3123

L-7200  L-7250

### SYRINGE (0.1 mL)



Dosierpritze (Nur Spritze)  
Bestell Nr. 080-1591

L-7200  L-7250

### PLUNGER TIP (0.1 mL)



Spritzendichtung  
Bestell Nr. 080-3280

L-7200  L-7250

### TEFLON LOOP (0.1 mL)



Probenschleife  
Bestell Nr. 810-3318

L-7200  L-7250

### SYRINGE (0.5 mL)



Dosierpritze (Nur Spritze)  
Bestell Nr. 080-1064

L-7200  L-7250

### PLUNGER TIP (0.5 mL)



Spritzendichtung  
Bestell Nr. 080-3276  
5 Pro Pack 080-3387

L-7200  L-7250

### TEFLON LOOP (0.5 mL)



Probenschleife  
Bestell Nr. 810-3319

L-7200  L-7250

### SYRINGE (1.0 mL)



Dosierpritze (Nur Spritze)  
Bestell Nr. 080-1592

L-7200  L-7250

### PLUNGER TIP (1.0 mL)



Spritzendichtung  
Bestell Nr. 080-3277

L-7200  L-7250

### TEFLON LOOP (1.0 mL)



Probenschleife  
Bestell Nr. 810-3320

L-7200  L-7250

### SYRINGE (2.5 mL)



Dosierpritze (Nur Spritze)  
Bestell Nr. 080-1092

L-7200  L-7250

### PLUNGER TIP (2.5 mL)



Spritzendichtung  
Bestell Nr. 080-3278

L-7200  L-7250

### TEFLON LOOP (2.5 mL)



Probenschleife  
Bestell Nr. 810-3321

L-7200  L-7250

# Autosampler Model LaChrom Classic L-7200 & L-7250

## SYRINGE (5.0 mL)



Dosierpritze (Nur Spritze)  
Bestell Nr. 080-1593

L-7200  L-7250

## PLUNGER TIP (5.0 mL)



Spritzendichtung  
Bestell Nr. 080-3279

L-7200  L-7250

## TEFLON LOOP (5.0 mL)



Probenschleife  
Bestell Nr. 810-3322

L-7200  L-7250

## NEEDLE (Standard)



Nadel Standard  
Bestell Nr. 810-3030

L-7200  L-7250

## NEEDLE GROOVED



Nadel mit Nut  
Bestell Nr. 810-3329

L-7200  L-7250

## NEEDLE SHARP EDGED TAP



Nadel angeschrägt  
Bestell Nr. 810-3110

L-7200  L-7250

## NOZZLE JOINT



Probennadel Halterung  
Bestell Nr. 080-1817

L-7200  L-7250

## NEEDLE SLEEVE



Ferrule Probennadel  
Bestell Nr. 635-6247

L-7200  L-7250

## TUBE ASSY



Schlauch Spülflüssigkeit  
Bestell Nr. 810-3090

L-7200  L-7250